

Flexometer

DOLIs Compression Flexometer

Since 1985 DOLI produces the fully automated Flexometer know also as "GOODRICH FLEXOMETER".

First of all, the automatic sample handling was important to our customers, followed by the height measurement and cooling stations, an automatic temperature chamber incl. preheating station and a sample silo for 20 samples.

Depending on the construction type of the Flexometer described in ASTM, ISO, or DIN, the frequency has been restricted to 30Hz, the load is defined by weights and the stroke is changed manually.

Of course DOLI continues developing the Flexometer.
See our overview!

Here you will find your overview folder with further information.